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- TriAuto ZX
- Tri Auto ZX2
- DentaPort
- TriAuto mini

EN

## Preparation of transfer instruments (turbines, scalers, straight and contra angle handpieces)

Practice hygiene has a key role to play in dentistry. As a responsible dental company, we assist dental surgeons in complying with directives and legislation established within the scope of dental hygiene.

The guidelines on the following pages guide you through the details with regard to the process of preparing your turbines, air motors, air scalers, straight handpieces and contra angles for micro motors and endodontic contra angle handpieces. We accompany you from the completion of treatment via cleaning to reuse.

But it is also important to us that you gain long-lasting pleasure from the high-quality products from Morita. This requires regular care, especially the high-performance turbines from the TwinPower series. TwinPower turbines are ceramically embedded. This is why it is absolutely necessary to oil the turbines with Morita Multi Spray at least twice daily or after a maximum of 20 minutes of operation.

Depending on the intensity of use, the turbines must be oiled more frequently if necessary. This interval arises automatically within the scope of routine re-preparation after every patient. If the turbines are only oiled irregularly or not at all, this leads to permanent damage and premature breakdown of the turbine. This results in consequential costs in every case.


Please also read the extensive instructions for use of the individual instruments. Pay attention to the warnings and tips in the instructions in particular.


## 1. External cleaning and disinfecting

Perform directly after treatment


### Take protective measures for possible infections

- Always wear protective gloves and mouth and eye protection during preparation.
- Take into consideration national regulations and directives (e.g. ISO 17665, UVV- VBG (Accident Prevention Regulation — Regulations of the Professional Associations) 103, RKI, etc.)
- Pay attention to the manufacturer's instructions on cleaning and disinfection products.


### Removing grinding instrument

- Remove grinding instrument (drill) from the collet. Thumb forceps are suitable for this to prevent injury.
- Place grinding instrument (drill) into a suitable disinfecting solution. Pay attention to the manufacturer's instructions on care of the grinding instrument.

### Remove handpiece from the tube


- Withdraw turbine from the quick coupling.
- Remove straight / contra angle handpiece from the ISO coupling on the motor.

### Cleaning

- Clean turbine or straight and contra angle handpiece under hot running water.
- It is best to use a soft brush.

Do not place in baths of disinfectant and do not clean in an ultrasonic bath.


### External disinfecting

- Thoroughly wipe turbine or straight and contra angle handpiece with a fiber-free disinfecting cloth. An appropriate disinfecting spray can be used instead of the fiber-free disinfecting cloth.

## 2. Internal cleaning and disinfecting

Perform within 30 minutes of last use. Choose the manual method 2.1 or alternatively the mechanical method 2.2.

Morita turbines, straight and contra angle handpieces with this symbol are suitable for both mechanical cleaning and disinfecting with a thermal disinfectant (RDG), and for manual internal cleaning and disinfecting with WL-clean and WL-cid.

## 2.1 Manual internal cleaning and disinfecting with WL-clean and WL-cid


(Manufacturer ALPRO MEDICAL) Pay attention to the manufacturer's instructions on care products

### Non-protein fixing cleaning with WL-clean

- Fit a suitable adapter to the WL-clean tin.
- Then place the turbine or straight and contra angle handpiece from the drive side onto the adapter.
- Spray through the turbine or straight and contra angle handpiece using three puffs of WL-clean lasting 2 seconds each.
- Remove the turbine or straight and contra angle handpiece from the adapter and set aside.
- Allow the cleaning agent to act for at least 1 minute.

## 2.1 Manual internal cleaning and disinfecting with WL-clean and WL-cid

### Disinfecting with WL-cid

- Fit a suitable adapter to the WL-cid tin.
- Then place the turbine or straight and contra angle handpiece from the drive side onto the adapter.
- Spray through the turbine or straight and contra angle handpiece with WL-cid for approximately 3 to 5 seconds.
- Remove the turbine or straight and contra angle handpiece from the adapter and set aside.
- Leave the disinfectant to act for approximately 1 to 2 minutes.

### Blowing out with WL-dry / WL-Blow

- Place a suitable adapter onto the WL-dry tin or connect the WL-Blow to the appropriate quick coupling on the compressed air hose.
- Place the turbine or straight and contra angle handpiece from the drive side onto the adapter.
- Blow compressed air through the turbine or straight and contra angle handpiece for approximately 3 seconds.
- Remove the turbine or straight and contra angle handpiece from the adapter.

## 2.2 Mechanical internal cleaning and disinfecting with a thermal disinfecter (RDG)

RDG can be used instead of manual internal cleaning. Place the assembled transfer instruments on the relevant adapters or in a dismantled condition into the RDG. Ensure that a suitable instrument adapter is used and the filter plates have been replaced at the prescribed intervals.

Pay attention to the manufacturer's instructions on the RDG devices.

## 3. Care with oil spray

At least on a daily basis at midday and in the evenings as well as before any sterilization and after any thermal disinfecting. Oiling the ball bearing is vital to ensure the durability of the instrument.

The TwinPower high-performance turbine in particular requires careful and regular care with the appropriate lubricants. Choose the manual method 3.1 or alternatively the mechanical method 3.2.


### **Morita Multi Spray**

Thanks to its high level of cohesion, Morita Multi Spray oil is ideally suited to pneumatic instruments (turbines, air scalers, pneumatic motors). It can cope with the high air throughput associated with drive air and is not „blown out“ from ball bearings.

### 3.1 Manual care with the Morita spray adapter

Morita Multi Spray oil is also highly fluid, which makes it very suitable for maintenance of low-speed instruments such as ISO straight and contra-angle handpieces. It will flow in all corners of the transmission. Additives ensure that the oil film does not break even in heavy use. The good lubricating properties (light running) ensure that less power is lost in transmission and is therefore not converted into useless heat. The handpieces are not as fast hot and almost the total torque of the micromotor comes to the bur. This means Morita Multi Spray is perfectly suitable for any heavy-duty rotating instruments you may use.

#### Oiling the ball bearing

- Shake oil spray can before use.
- Fit a suitable adapter onto the tin.
- Press down nozzle lever.
- Insert nozzle into the turbine or straight and contra angle handpiece from the drive side.
- Spray the instrument with oil spray for 2 to 3 seconds.
- Remove the instrument from the adapter.

**Never direct the spray nozzle at a person and never use the spray can in proximity to a naked flame.**

### 3.1 Manual care with the Morita spray adapter

#### Oil the collet at least once a week

- Lift the lever on the nozzle connection and swivel the chuck nozzle into position.
- Cover the instrument head with gauze. This protects the surrounding area from oil pollution.
- Insert the nozzle into the chuck and spray for 2 to 3 seconds.


### 3.2 Mechanical care with the Lubrina care device

As an alternative to manual care, a care device (e.g. Lubrina) can be used. Pay attention to the manufacturer's instructions on care devices.


Lubrina by Morita has an integrated device for collet care. For competitors' care devices without a mechanism for oiling the collet, this must be performed manually (see above).

### 4. Check the FG collet


#### Check the FG collet once a week according to RKI

- Insert the test instrument into the drill chuck.
- Clamp the pin as for a grinding instrument (drill).
- Grip the instrument head with your hand. Do not touch the instrument handle piece. This could result in damage to the collet.
- Pull apart the axial instrument head and test instrument.

#### Blue line visible

- The necessary holding force is attained when the blue ring is visible on the shank. If the FG shank disengages after the blue ring, the collet is fine and the instrument can continue to be used.

#### Blue line not visible

- If the FG shank disengages before the blue ring is visible, the collet no longer has sufficient tension. If this occurs, the cartridge must be replaced.

Instruments which have not been used surgically and are also not intended for this use can now be used again or stored in a dust-protected environment (continued in 6.).

- Semi-critical medical device: medicinal products coming into contact with mucous membranes or pathologically altered skin.
- Medical device semi-critical B: e.g. transfer instruments for general, restorative or orthodontic treatment.

## 5. Sterilization:

Instruments intended for use in the surgical field are shrink-wrapped into appropriate sterile packaging and sterilized in B-class autoclaves (EN DIN 13060).

- Critical medical device: medical devices for the application of sterile pharmaceutical and medical devices which penetrate through the skin or mucous membranes and thus come into contact with blood, internal tissues or organs, including wounds.
- Medical device critical B: e.g. transfer instruments for surgical, periodontological or endodontic treatment.

## 6. Storage


### Sterilization in the autoclave

- Wipe the turbine or straight and contra angle handpiece thoroughly with a piece of gauze and then autoclave.
- The recommendations of the RKI from 4/2006 apply, in addition to the usage instructions of the manufacturer of sterile films and autoclaves.
- Steam sterilizers are permitted which conform to EN 13060 class B and which are authorized for the sterilization of straight and contra angle handpieces and turbines.
- Do not exceed a temperature of 135°.
- The instruments can be destroyed if they are not properly oiled before autoclaving.
- Do not leave instruments in the autoclave overnight.


- The instruments are stored, according to necessity, in low bacteria or sterile conditions, unpacked or packed, shielded from dust, clean and dry.
- Storage should take place observing the expiry dates for sterile reuse in accordance with the specifications from the sterile packaging manufacturer. Pay attention to the quality of the welded seam in particular.

### Blowing out

- Connect the turbine using the quick coupling on the instrument tube, the straight and contra angle handpiece onto the ISO coupling of the micro motor.
- Leave the instrument to run for approximately 15 seconds to blow out the superfluous oil.
- Wipe head with disinfecting cloth.

**Do not touch the push button when the instrument is operated to remove superfluous spray. Technical defects can arise in the case of non-compliance.**

### Morita chuck tester

Morita collet testing instrument (green dynamometer) for turbines and straight and contra angle handpieces:

Order no. N5010695


### Morita care spray

Morita Multi Spray for all turbines, contra-angles and also suitable for endodontic contra angle handpieces.

Order no. 2698

## 9. Transfer instruments

**Morita spray adapter**

Morita spray adapter for straight and contra angle handpieces:  
 Morita turbine couplings:  
 Sirona turbine couplings:  
 KaVo turbine couplings:  
 NSK turbine couplings:  
 W+H turbine couplings:

Order no. N5071345  
 Order no. N5010307  
 Order no. N5010297  
 Order no. N5010293  
 Order no. N5010299  
 Order no. N5010295

Morita instruments with this symbol are suitable for both mechanical cleaning and disinfecting with a thermal disinfectant (RDG), and for manual internal cleaning and disinfecting with WL-clean and WL-cid.

**Pneumatically operated instruments**

- All TwinPower turbines
- AirTorx air motor
- AirSolfy air scaler

**Handpieces and contra angle handpieces for micromotors**

- All straight and contra angle handpieces in the Tokyo series
- All straight and contra angle handpieces in the TorqTech series

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