# Reprocessing Beach Cannulae

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EN

Reprocessing Beach Cannulae  
(Re-usable suction cannulae from Morita)

## Reprocessing Beach Cannulae

(Re-usable suction cannulae from Morita)

Practice hygiene plays a major role among dentists. As a dental company which is aware of its responsibility, we help practitioners adhere to the directives and legislation they face in respect of dental practice hygiene.

The guidelines on the following pages will lead you through the details on the preparation procedures for re-using your Beach suction cannulae. We will assist you all the way from the completion of treatment to cleaning and re-using the cannulae.

## Preparation instructions

Beach suction cannulae are available both as disposable products made of plastic as well as re-usable metal cannulae. Re-use is only possible with metal cannulae and their blue rubber tips.

The preparation of medical products for safe re-use is on the assumption that they fall under one of the risk groups that are shown on page 30. Mechanical and/or manual treatment processes are available depending on the type of product and risk evaluation.

Medical product	Type of application
<b>Non-critical</b>	Only contact with intact skin.
<b>Semi-critical A:</b> without particular requirements on treatment <b>Semi-critical B:</b> with particular requirements on treatment	Contact with the mucosa or abnormally modified skin.
<b>Critical A:</b> without particular requirements on treatment  <b>Critical B:</b> with particular requirements on treatment	Penetration of the skin and/or mucosa and at the same time contact with blood, internal tissue or organs.  Contact with wounds. Medical products for use with blood, blood products and other sterile medications.

## Preparation instructions


Please observe the respective national regulations for the cleaning, disinfection, sterilization and storage of Beach cannulae.

Please also read the extensive instructions for use of the individual instruments. Pay particular attention to any warnings and comments they contain.


## 1. External cleaning and disinfection

Carry out immediately prior to treatment.


### Take appropriate protective measures against possible infection

- Always wear protective gloves as well as mouth and eye protectors during preparation treatment.
- Observe national regulations and directives (e.g. ISO 17665, UVV-VBG 103, RKI, ...).
- Observe manufacturer's instructions for cleaning and disinfection products.


### Disconnect suction cannula(e) from handpiece

- Remove larger amounts of contamination with a disinfecting one-way cloth.
- Pull the spray mist cannula and rubber tip off the large suction handpiece separately.
- Pull the saliva ejector cannula off the small saliva ejector handpiece.


### Transport/Storage

- After use, collect Beach cannulae and rubber tips in a suitable, enclosed container filled with instrument disinfection solution and transport safely to disinfection point.
- Observe manufacturer's instructions for cleaning and disinfection products. Rubber tips suffer in particular if the disinfection solution is mixed incorrectly.

### Cleaning

- Clean the Beach cannulae and rubber tips under flowing warm water.
- For best results use a soft brush.

### External disinfection

- Wipe off the cannulae and rubber tips completely using a lint-free disinfection cloth. A suitable disinfection spray may be used as an alternative to lint-free cloth.

## 2. Internal cleaning and disinfection

Carry out within 30 minutes of last use. Choose either the manual method 2.1 or, alternatively, the mechanical method 2.2.

Beach suction cannulae (Morita) made of metal and rubber tips are suitable for both mechanical cleaning and disinfection with a thermo-disinfector, as well as manual internal cleaning and disinfection.

## 2.1 Manual cleaning and disinfection by wet-chemical process

- Observe care product manufacturer's instructions.
- Use VAH/DGHM-listed instrument disinfection and cleaning solutions.
- Insert suction cannulae and rubber tips completely.
- Check the cannulae and rubber tips for residual contamination upon completion of the cleaning and disinfection cycle. Repeat the cycle if necessary.
- Rinse cannulae and rubber tips with water.
- Dry with disposable cloths or using compressed air.

## 2.2 Mechanical cleaning and disinfection with a thermo-disinfector

A thermo-disinfector unit may be used as an alternative to manual internal cleaning.

- Insert the cannulae in the thermo-disinfector on appropriate adapters.
- Pay attention that an appropriate instrument adapter is used and that the filter panels are replaced in regular cycles in compliance with regulations.
- Observe the manufacturer's instructions for thermo-disinfection equipment.
- Check the cannulae and rubber tips for residual contamination upon completion of the cleaning and disinfection cycle. Repeat the cycle if necessary.

### 3. Sterilization

Instruments that are intended for use in surgical areas are welded in suitable sterile packaging and sterilized in B-class autoclaves (EN DIN 13060).

- Critical medical products: Medical instruments for application of sterile medications and medical products that penetrate the skin or mucosa and at the same time come into contact with blood, internal tissue or organs, including wounds.
- Medical products, critical B: e.g. transmission instruments for surgical, periodontal or endodontic treatment.


**Instruments that have not been used for surgical purposes and are not intended for such use either, can now be re-used again or stored and protected against dust (go to Section 4.).**

- Medical products, semi-critical, that come into contact with mucosa or abnormally modified skin.
- Medical products, semi-critical B: e.g. transmission instruments for general, restorative or orthodontic treatment.

### Sterilization in autoclaves

- Wipe off suction cannulae and rubber tips completely with a piece of gauze and then autoclave.
- Observe the recommendations of the Robert Koch-Institut of 4/2006 and the sterile foil and autoclave manufacturers' operating instructions.
- Steam sterilizers are permitted that comply with EN 13060 Class B and are approved for the sterilization of dental instruments.
- A temperature of 135° may not be exceeded.
- Observe process parameters.
- Avoid metal contact of sterile products with other instruments or the edge of the tray.
- Observe manufacturer's instructions for the correct placement of sterile products.
- Do not leave instruments in the autoclave overnight.

Observe all national and local directives and principles for the sterilization of medical products.



## 4. Inspection and functional testing

Use and frequent reprocessing has effects on the condition of both the cannulae and the rubber tips.


- The end of a product's service life is determined by wear and visible damage as well as color changes.
- Check cannulae and rubber tips and replace if necessary.
- Disinfect and sterilize damaged cannulae and rubber tips before disposing.


## 5. Storage


- Depending on requirements, store suction cannulae and rubber tips in low-germ environments or sterile, packed or unpacked, dust-protected, clean and dry.
- Store under adherence to the use-by dates for sterile re-use according to the manufacturer's instructions for the sterile packaging. Particular attention is to be paid to the quality of the welded seam.

- Attach the suction cannulae and rubber tips to the suction handpieces of the treatment unit immediately prior to commencement of patient treatment.

Order no.	Form	Description
2087		Rubber tip "B" for saliva ejector cannulae: <ul style="list-style-type: none"> <li>■ 2079</li> <li>■ 2088</li> </ul>
2086		Rubber tip "A" for cannulae: <ul style="list-style-type: none"> <li>■ 2080</li> <li>■ 2081</li> <li>■ 2082</li> <li>■ 2083</li> </ul>


Order no.	Form	Description
2090		Saliva ejector cannulae "C" for saliva ejector hose
2079		Suction cannulae for saliva ejector hose
2088		Saliva ejector cannulae "C-N" ■ adaptable direct on the suction handpiece

Order no.	Form	Description
2084		Suction cannulae "B" ■ surgical ■ suitable for cannulae 2081 or 2083
2089		Suction cannulae "B-N" ■ surgical ■ adaptable direct on the suction handpiece
2080		Suction cannulae "A1-120" ■ 120° curved ■ with mesh screen

Order no.	Form	Description
2081		Suction cannulae "A2-120" <ul style="list-style-type: none"> <li>■ 120° curved</li> <li>■ without mesh screen</li> </ul>
2082		Suction cannulae "A1-90" <ul style="list-style-type: none"> <li>■ 90° curved</li> <li>■ with mesh screen</li> </ul>

Order no.	Form	Description
2083		Suction cannulae "A2-90" <ul style="list-style-type: none"> <li>■ 90° curved</li> <li>■ without mesh screen</li> </ul>

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