

■ Variation

Tray arm	Over-the-patient tray
Dental chair	Contour type
Seat leather color variations	G10 Shiny Silver (G10 SS) G10 Splash Yellow (G10 SY) G10 Honey Pink (G10 HP) G10 Fresh Green (G10 FG) G10 Sweet Brown (G10 SB) G10 Cool Violet (G10 CV) G10 Gloss Black (G10 GB)

■ Specifications

Name	SIGNO G10
Model	TU101
Input voltage	Adjusted to local line voltage
Power consumption	1.5kVA
Patient seat elevation system	Hydraulic motor
Seat height	750mm MAX 400mm MIN
Backrest angle (from horizontal line)	Raised 70° Lowered 2°
Weight	Approx. 240kg

- The product specifications and appearance are subject to change without notice.
- The printed product colors may differ from the actual colors of the product.
- The photographs in this catalog show the dental chair with some optional devices.



J.MORITA CORPORATION
3-33-18 Tarumi-cho Suita-City, Osaka, 564-8650 Japan
Tel:+81-6-6380-1521

J.MORITA TOKYO MFG.CORP.
7129 Komuro, Ina-Machi, Kitaadachi-gun, Saitama-Ken, 362-0806 Japan
Tel:+81-48-723-2621



Signo G10

Evolution of a new era

The Signo G10 is an evolution of the proven Signo range of dental treatment units from Morita. Signo G10 offers a contemporary streamlined aesthetic appeal, complemented by a range of advanced dental elements designed to meet the requirements of the modern dental professional.


LED Light



LED Operating Light

Luna Vue ES – Signo G10 offers a new integrated LED operating light which provides vivid natural color with high luminescence, while still providing practical advantages of high energy efficiency and long lasting durability.

Seat



Comfortable Patient Seating

The design of the Signo G10 backrest and seat ensures uncompromising patient comfort and support. Designed to ensure a feeling of security for the patient, the Signo G10 delivers maximum comfort. The leather upholstery is available in a range of colours, all incorporating a modern cross-hatched design texture ensuring esthetic appeal, as well as comfort.

Instrument Holder



Instrument Holder and Tray

The Signo G10 instrument holder provides a seamless transition to the generously sized operator's tray. The instrument holder incorporates a flexible array of instruments including a three-way syringe, highspeed handpieces, an electric micromotor and scaler. Optionally available with TwinPower or TorqTech handpieces.

Basin



Glass basin

Focused on quality and cleanliness the Signo G10 offers a contemporary glass basin design. The basin is easily removed for cleaning and is available in a range of colours to complement your upholstery choice.

Tray Arm



Tray Arm

The operator's tray and instrument holder support arm is designed to provide maximum flexibility to allow easy patient access to the treatment area. The operator can easily move the tray vertically and horizontally to match their desired operating position. While an arm movement lock ensures a safe environment for both patient and operator.

Chair



Chair design

In addition to contemporary aesthetic appeal, the Signo G10 chair has been designed to provide maximum flexibility for patient access. The height of the chair can be adjusted from a minimum of 400mm to a maximum of 750mm. This wide range ensures easy access for patients of all ages.

Basin Color



Clear



Green



Orange



Pink

Leather Color



G20 Gloss Black (G20 GB)



G20 Shiny Silver (G20 SS)



G20 Fresh Green (G20 FG)



G20 Cool Violet (G20 CV)



G20 Sweet Brown (G20 SB)



G20 Splash Yellow (G20 SY)




G20 Honey Pink (G20 HP)

Standard



Three-way syringe
(WS-201)



Vacuum syringe VS14WD2
(with suction power control valve)




Luna Vue ES Light
(Single LED)


Option




Highspeed handpiece
TwinPower (with optic)
TwinPower Ultra Series (with optic)



Electric micromotor
TORX TR-91, TR-91-O (with optic)



AIR TORX TRA-200
Straight handpiece(ST21U)
Contra shank(US121L)
Contra head(QL26B)



Saliva ejector
SE14S
(with suction power control valve)



Assistant tray




Luna Vue EL Light
(12 LED)



EL mirror



Headrest
(Articulated headrest)



Movable armrest
At the dentist's and assistant's side



Movable armrest
At the assistant's side



Small tray mounted on light pole
Tissue box holder
(Over arm type)



Bacteria filters
(water and air)



Foot controller
<Extended function with lever>
(FC201A)



Triple viewer (PFV95)



Operating stool LN (OS-LN)
Width 510mm(seat width 400mm)
Depth 540mm(seat width 450mm)
Seat Height 445~635mm
Height Adjustable Range 190mm
Weight Approx.16kg